



Model:
Phase Contrast set

Description:
Phase contrast set with IOS W-PLAN PH objective 40x

Condenser (M-1124)	1.25 N.A. Fixed frontal lens. Iris aperture diaphragm, plus slider with phase contrast annuli. Positions: Brightfield, 10x phase contrast, 40x phase contrast.
Objective (M-1122.N)	Negative Phase Contrast 40x/0.65 N.A. Planachromatic IOS (infinity-corrected tube length). Working Distance: 0.54 mm (spring-loaded). Cover Glass Thickness: 0.17 mm.
Accessories	<ul style="list-style-type: none"> - Centering Telescope (M-1004.N) - Green filter (M-977) - Hard case