



Model:

MS-2

Typology:

EDUCATIONAL STEREOMICROSCOPE

Description:

Entry- level stereomicroscope for educational field.

Head	Monocular, 45° inclined, 360° rotating.
Eyepieces	Wide Field, WF 10x/16.
Objective	2x fixed. Working distance: 57 mm.
Stand	Pillar with focus mechanism.
Focusing	By means of a rack and pinion controlled by a pair of knobs placed on both sides of the microscope.
Stage	Fitted with specimen clamps.
Illumination	Without illumination (external illuminators are available).
Dimensions	HEIGHT: 255 mm WIDTH: 108 mm DEPTH: 147 mm WEIGHT: 1.1 kg
Accessories included	Manual and dust cover.