



*Large stand, pillar type, with transmitted LED illumination*

<b>Focusing System</b>	Type	Coarse
	Coarse total travel (mm)	50
	Adjustable tension	Yes
	Head holder internal diameter (mm)	76
<b>Stand</b>	Type	Pillar
	Pillar diameter (mm)	32
<b>Stage</b>	White diffusion plate for transmitted ill.	Yes, diameter 82 mm (diameter for accessories holder 94 mm)
<b>Transmitted Illumination</b>	Type	LED
	Light source power (W)	2
	Illuminance (lux)	22,000
	Brightness control	Manual
	Lifetime (hours)	> 65,000
	Temperature (K)	6,300
	Max. required power (W)	3
<b>Power Supply for Illumination</b>	Type	Internal
	Microscope connector	C13 plug
	Power plug type	Schuko
	Input voltage	100/240 Vac, 50/60 Hz
<b>Accessories Included</b>	Dust cover	Yes
	User Manual	Digital version (downloadable)
<b>Product Dimensions</b>	Total height (mm)	410
	Base height (mm)	50
	Base width (mm)	320
	Base depth (mm)	290
<b>Product Weight</b>	(kg)	4.2