



Advanced overhanging stand

<b>Stand</b>	Type	Hinged overhanging
	Pillar diameter (mm)	32
<b>Accessories Included</b>	Dust cover	Yes
	User Manual	Digital version (downloadable)
<b>Additional Information</b>	Maximum sample height (mm):	440 (from bench)
<b>Product Dimensions</b>	Total height (mm)	425
	Base height (mm)	50
	Base width (mm)	210
	Base depth (mm)	260
	Horizontal arm (mm)	515 (overall)
<b>Product Weight</b>	(kg)	18.8

