



Ball-bearing overhanging stand, with focusing system and head holder

Focusing System	Type	Coarse
	Coarse total travel (mm)	50
	Adjustable tension	Yes
	Head holder internal diameter (mm)	76
Stand	Type	Simple overhanging
	Pillar diameter (mm)	32
Accessories Included	Dust cover	Yes
	Allen wrench	Yes
	User Manual	Digital version (downloadable)
Additional Information	Maximum sample height (mm):	270 (from bench)
Product Dimensions	Total height (mm)	430
	Base height (mm)	58
	Base width (mm)	210
	Base depth (mm)	255
	Horizontal arm (mm)	790
Product Weight	(kg)	16.3

