



*Hinged overhanging stand, with focusing system and head holder*

<b>Focusing System</b>	Type	Coarse
	Coarse total travel (mm)	50
	Adjustable tension	Yes
	Head holder internal diameter (mm)	76
<b>Stand</b>	Type	Hinged overhanging
	Pillar diameter (mm)	32
<b>Accessories Included</b>	Dust cover	Yes
	Allen wrench	Yes
	User Manual	Digital version (downloadable)
<b>Additional Information</b>	Maximum sample height (mm):	440 (from bench)
<b>Product Dimensions</b>	Total height (mm)	430
	Base height (mm)	58
	Base width (mm)	210
	Base depth (mm)	255
	Horizontal arm (mm)	670
<b>Product Weight</b>	(kg)	19.8