



*Pantograph stand with head holder and focusing mechanism. Provided with table clamp*

<b>Focusing System</b>	Type	Coarse
	Coarse total travel (mm)	50
	Adjustable tension	Yes
	Head holder internal diameter (mm)	76
<b>Stand</b>	Type	Pantograph, with table clamps
	Pillar diameter (mm)	35
<b>Accessories Included</b>	Dust cover	Yes
		Yes
	User Manual	Digital version (downloadable)
<b>Product Dimensions</b>	Total height (mm)	330
	Total depth (mm)	880
	Horizontal arm (mm)	820
<b>Product Weight</b>	(kg)	4,8

