



Model:
SZN-T

Typology:
ROUTINE STEREOZOOM HEAD

Description:
Stereomicroscope head with zoom for routine and industrial applications.

Illumination	<p>The microscope has no illumination system, but it can be used with external illuminators (halogen, LED or optical fibers) or with the following models:</p> <ul style="list-style-type: none"> • ST-151 • ST-153 • ST-155 • ST-156 • SZ-STL8 • SZ-STL5LED
Focusing	<p>The microscope can be used with the following models of stand:</p> <ul style="list-style-type: none"> • ST-150 • ST-151 • ST-152 • ST-153 • SZ-STL1 • SZ-STL1H • SZ-STL2 • SZ-STL2H • ST-155 • ST-156 • SZ-STL8 • SZ-STL5 • SZ-STL5LED • SZ-STLX <p>As an alternative it's possible to use focusing systems with a 76 mm diameter head holder.</p>
Head	<p>Trinocular 360° rotating and 45° inclined. Adjustment +- 5 diopters on both eyepieces. Interpupillary adjustment 55-75 mm.</p>
Eyepieces	<p>Extra Wide Field, EWF 10x/22.</p>
Objective	<p>Parfocal achromatic objective 1X Working distance: 110 mm Zoom range: 0.67x - 4.5x (zoom factor: 6,71:1). Maximum magnification in standard configuration: 45x (can be increased to 225x).</p>
Accessories	<p>Instruction manual and dust cover included.</p>